Product Specifications

Name: MEJES IX E-ion removes Speckle Cosmetic Device
Model: MEJES9701FR00913
Charging Type: USB Charging
Power Voltage: 5V⎓1A
Power Output: 5W
Battery Capacity: 1800mAh
Charging Duration: 4 hours
Battery Working Time: 90 Minutes
Standby Time: 90 Days (If the unit hasn't been used for more than 90 days, please make it fully charged before use)
Size: 163×43×41mm

Product Structures

① Massager
② Gear
③ Power
④ Charge

Use

Press and hold for 2 seconds to turn on/off.

Press the ⇳ key to adjust the gear.

Note: It is recommended to use the lowest level for sensitive skin.

1.Clean the face and keep face moist.

2.Long press power to start the device.

3.Moisten the skin with negative E-Water

The massage panel gently touches the skin

Slowly moving on the skin

Important:

1. Keep skin moist with conductive gel .
2. After use , cleaning is not allowed within 4 hours
3.Negative E-Water Usage order: 20c-100c-500c-1000c... gradually

Safety Warnings

- Please read all instructions carefully for safety and proper use.
- Use only as directed.
- Immediately stop using if any malfunction or unusual reaction to treatment occurs.
- Precautions listed in this manual are very important. For your safety, please follow them carefully.
- Precautions are classified as WARNING or CAUTION in order to clearly specify the level of damage and possible injury or damage that may occur if the device is misused.

WARNING & CAUTION	Improper use of this device may result in death or injury
Prohibited	Do not disassemble, attempt to repair or modify this device in any way.
Remove Power Cord	Hold the head of the power cord instead of any part of the cord when unplugging the power cord from USB interface. Damage to the cord may cause electrical shock, short circuit or fire hazard.
	Always unplug the power cord from USB interface when not use. Do not handle the AC adapter with wet hands. Which may result in electrical shock, short circuit or fire hazard.
	If any malfunction occurs, unplug the power cord from the device immediately. Possible malfunctions may include button failure, smoke, unusual noise or smell, please stop using it.
Prohibited	Do not use if the power cord is damaged or if the wall socket is loose which may cause electrical shock, short circuit.
	Disabled for children under 17 years of age,Pregnant women are disabled.If you are sick, please follow the doctor's instructions.
	Do not use extension cords. This may cause electrical shock, overheating or fire hazard.
Prohibited	Not for use by children, elderly, or disabled person. Improper or unattented use may cause injury.

Warranty

This product is warranted for 1 year from the purchase date.
1.In accordance with the instructions to use, in the warranty period, it it's not work, factory will offer free maintenance. (Please submit the goods repair and warranty card)
2.Even in the warranty period, if it's damaged caused by the following circumstances,.
2.1 Malfunctions and damage caused by incorrect use
2.2Due to their own improper repair or transformation of the failure and damage
2.3After the purchase due to transport, drop, etc. arising from the failure and damage
2.4Occurrence of fire, earthquake, flood, thunder and other weather variations, abnormal voltage generated by the failure and damage
3.The warranty card will not be valid once it's separately from the user manual, please keep it safely.

Product Warranty Card

Order Number			
Purchase Name		Phone Number	
Purchase Date		Email	
Purchase Address			

Maintenance Record			
Order Number	Maintenance Part	Maintenance Result	Maintenance Person

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